



PRODUCT INFORMATION

Rapid Animal Total RNA Extraction Kit

Product information for AT4181 & AT4182:

Components	AT4181 50 Preps	AT4182 250 Preps
Solution A	50 ml	250 ml
Protocol	1	1

Introduction:

This kit uses a reagent similar to TRIzol for rapid RNA isolation from animal tissue, cells, and some of the plant tissue.

Feature:

- ✓ Rapid and handy: the whole procedure takes approximately 10 minutes.
- ✓ Very low genomic DNA contamination (no detectable by RT-PCR).
- ✓ Flexible, ideal for most animal tissue and cell culture.

NOTE: Care must be taken when working with RNA. It is important to maintain an RNase-free environment starting with RNA sample preparation and continue through purification and analysis. Use RNase free tubes, tips, gels. Wear gloves at all times.

Procedure:

1. Preparation
 - A. Estimate the amount of Solution A to be used.



For 1 ml of Solution A, the maximum amount of tissue to be used is 100 mg; the maximum amount of cells to be used is 1×10^7 .

- B. Prepare enough chloroform, isopropanol, 75% ethanol and RNase-free water.
2. Homogenization
 - A. *For adherent cells*
Discard the culture medium, add 1 ml of Solution A to ~10 cm² cultured cell, and mix gently by pipetting.
 - B. *For suspension cells*
Collect cells by centrifugation, discard the supernatant and add 1 ml of Solution A to $1-5 \times 10^6$ cells, and mix gently by pipetting. The amount of cells should not exceed 1×10^6 for fibroblasts or carcinoma cell.
 - C. *For fresh animal tissue*
Cut the tissue into pieces; add 1 ml of Solution A to 50-100 mg tissue, homogenate for 30 seconds using a homogenizer.
 - D. *For Samples in RNAlocker*
Remove the liquid on the surface prior to cut, and then treat as fresh tissue.
3. Fraction Separation
 - A. Add 0.2 volume of chloroform to the lysate, (eg. 0.2 ml chloroform for 1 ml lysate), vortex for 30 sec.
 - B. Centrifuge at no more than 12,000 x g (13,000 -15,000 rpm) for 3 minutes at room temperature.
 - C. Transfer the supernatant into a new 1.5 ml centrifuge tube.
 - D. Repeat step 3 for 1-2 times for tissues rich in lipid.
4. Precipitation
 - A. Add equal volume of isopropanol to the supernatant, vortex for 30 seconds.
 - B. Centrifuge at no more than 12,000 x g (13,000-15,000 rpm) for 3 minutes at room temperature, deposit may form at the bottom of the centrifuge tube.



Bio Basic Inc.

A world leader in serving science

-
- C. Discard the supernatant carefully.
 5. Primary Wash
 - A. Add 1 ml of 75% ethanol to the RNA deposit, vortex for 30 seconds.
 - B. Centrifuge at no more than 7,500 x g (13,000-15,000 rpm) for 3 minutes at room temperature.
 - C. Discard the supernatant carefully.
 6. Secondary Wash
 - A. Repeat step 5 once.
 - B. Short spin at no more than 7,500 x g (13,000 rpm), remove the residual ethanol carefully.
 7. Dissolving the RNA
 - A. Incubate the tube with lid open at room temperature for 1-2 minutes.
 - B. Add 50-100 µl of RNase-free water to dissolve the pellet, store RNA solution at -80°C.
 8. RNA integrity, yield and purity can be analyzed via electrophoresis and spectral analysis.

Storage:

Transport at room temperature, store all components at 4°C upon received.



Bio Basic Inc.

A world leader in serving science

**PRODUCTS ARE INTENDED FOR BASIC SCIENTIFIC
RESEARCH ONLY!
NOT INTENDED FOR HUMAN OR ANIMAL USE!**

Please visit www.biobasic.com



A World Leader in Serving Science