



A world leader in serving science

BIO BASIC INC.

BS4654 (4 preps)

BS466 (20 preps)

**Plasmid DNA Extraction
MAXI Prep Kit**

Version 3.0
ISO9001 Certified

20 Konrad Cres, Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com

Introduction

Bio Basic Plasmid DNA Extraction Maxiprep Kit is an excellent tool offering a rapid and economic method to purify plasmid DNA from bacterial cultures. This technology is based on alkaline lysis and purification by Anion-exchange chromatography. Compared with other harmful and time-consuming procedures such as phenol / chloroform extraction and ethanol precipitation, Bio Basic Plasmid DNA Extraction Kit shortens the handling time to about 2 hours. The high quality plasmid DNA can be used directly for any downstream application.

Specification

Sample Size	Yield	Handling Time
100-250mL of bacterial culture for high copy plasmids	Up to 500 µg for high copy-plasmids	About 2 hours
200-400 mL of bacterial culture for low copy plasmids		

Content

Component	BS4654 (4 preps)	BS466 (20 preps)
MAXI 1 <i>Resuspension Solution</i>	55mL	2X110 mL
MAXI 2 <i>Cell Lysis Solution</i>	55 mL	2X110 mL
MAXI 3 <i>Neutralization Solution</i>	55 mL	2X110 mL
MAXI 4 <i>Equilibration Solution</i>	54 mL	2X135 mL
MAXI 5 <i>Washing Solution</i>	4X55 mL	4x275 mL
MAXI 6 <i>Elution Solution</i>	54 mL	2X135 mL
RNase A (50 mg/mL)	88 µL	2X220 µL
MAXI Column	4 pcs	20 pcs
User Manual	1	1

Additional Materials Required



1. 50 mL centrifuge tube.
2. Isopropanol
3. 70% Ethanol

Important Notes

1. Solutions provided in this kit contain irritants, wear gloves and lab coat when handling.
2. Briefly spin RNase A tube to remove drops from the inside of the lid. Transfer tube contents into **MAXI 1 Resuspension Solution** bottle. Add 250 μ L of **MAXI 1 Resuspension Solution** into RNase A tube, rinse tube inside and transfer back into **MAXI 1 Resuspension Solution** bottle. Store at 4°C.
3. Check **MAXI 2 Cell Lysis Solution** before use. Warm **MAXI 2 Cell Lysis Solution** at 37°C if any precipitation formed. Prevent vigorous shaking of the **MAXI 2 Cell Lysis Solution**.
4. To avoid acidification of **MAXI 2 Cell Lysis Solution** from CO₂ in the air, close the bottle immediately after use.

Brief Procedure

Brief Procedure



General Protocol

1. Harvest the bacterial culture by centrifugation at 6,000 x g for 15 minutes.
2. Add 10 mL of **MAXI 1 Resuspension Solution** (RNase A added) and resuspend the cell pellet by vortexing or pipetting.
3. Add 10 mL of **MAXI 2 Cell Lysis Solution** and mix gently by inverting the tube 10 times. Do not vortex to avoid shearing of genomic DNA.
4. Incubate for 3 minutes at room temperature until lysate clears.
5. Add 10 mL of **MAXI 3 Neutralization Solution** and mix immediately by inverting the tube 10 times. Do not vortex to avoid shearing of genomic DNA.
6. Centrifuge at 15,000 x g for 20 minutes at 4°C.
7. Transfer supernatant from step 6 in a new tube. Centrifuge at 15,000 x g for 20 minutes at 4°C.
8. Place a **MAXI** Column into a 50 mL centrifuge tube, add 10 mL of **MAXI 4 Equilibration Solution** to equilibrate the **MAXI** Column and allow the column to empty by gravity flow. Discard the filtrate.
9. Transfer the supernatant from step 7 to the equilibrated **MAXI** Column, and allow column to empty by gravity flow. Discard the filtrate.
10. Add 25 mL of **MAXI 5 Washing Solution** to wash the **MAXI** Column and allow the column to empty by gravity flow. Discard the filtrate.

11. Repeat Step 10.
12. Place **MAXI** Column into a clean 50 mL centrifuge tube (not provided) and add 12 mL of **MAXI 6** Elution Solution to elute DNA by gravity flow.
13. Precipitate DNA by adding 9mL of isopropanol to the eluted DNA from Step 12.
14. Mix gently and centrifuge at 20,000 x g for 30 minutes at 4°C.
15. Carefully remove the supernatant and wash the DNA pellet with 5mL of room temperature 70% ethanol.
16. Centrifuge at 20,000 x g for 10 minutes at 4°C.
17. Carefully remove the supernatant and air-dry the DNA pellet for 10 minutes.
18. Dissolve the DNA pellet in a suitable volume of 10mM tris pH8.5 or ddH₂O.



Troubleshooting

Low yield

Bacterial cells were not lysed completely.

- Too many bacterial cells were used.
- After **MAXI 3 Neutralization Solution**, break up the precipitate by inverting.
- DNA pellet was lost after precipitation.
- DNA pellet was insufficiently redissolved.

Purified DNA doesn't perform well in downstream application

RNA contamination


- Make sure that RNase A has been added in **MAXI 1 Resuspension Solution** when first using. If RNase A added in **MAXI 1 Resuspension Solution** has expired, add additional RNase A.
- Too many bacterial cells were used, reduce the sample volume.
- Elution buffer contains EDTA.

Genomic DNA contamination

- Do not use overgrown bacteria culture.
- During **MAXI 2 Cell Lysis Solution** and **MAXI 3 Neutralization Solution** addition, mix gently to prevent genomic DNA shearing.
- Lysis time was too long (over 5 minutes).

Too much salt residual in DNA pellet

- Wash the DNA pellet twice with 70% ethanol.



**PRODUCTS ARE INTENDED FOR BASIC
SCIENTIFIC RESEARCH ONLY!
NOT INTENDED FOR HUMAN OR ANIMAL USE!**

Please visit www.biobasic.com



A world Leader in Serving Science